

The FG463IR-UGRO display upgrade kit is a cost effective way to upgrade any existing printer equipped **SR**® in-floor scale to include all the great new features that a new scale has to offer including auto weigh, zero tracking, auto print with auto cut printer and BMI functionality.



Touchless operation - The scale will automatically turn on, lock in a stable weight, and print out a record with weight, date, and time stamp. (BMI can also be calculated and added to the record)

Adjustable viewing angle - The display can rotate to adjust viewing angle. It can also be rotated to hide the display and buttons during normal operation or when not in use.

SYSTEM SPECIFICATIONS

Max Weight	1000 Lb /454 Kg	Auto Zero	Scale will zero track
Display Type	16 Character LCD w/LED Backlight	Auto Weigh/Print	The scale will automatically turn on, lock in a stable weight and print with auto cut
Display Resolution	0.1% +/-1 digit of displayed resolution	Memory	Last stable weight is automatically stored in memory. Press RECALL to retrieve
Display Units	Pounds / Kilograms	BMI Calculation	Press BMI and input Height. Press ENTER to view BMI information
Power	120VAC 50-60 Hz		

