

### Quality. Reliability. Accuracy.

# Gram Scales Model SR300/310/320/325

When your weighing needs require gram accuracy, the SR300 Series is the right tool for the job. Available in four different weight ranges, the SR300 Series can handle a variety of tasks from nutritional portion measurements to sub-gram measurements for tracking neo-natal diaper output.

Constructed of high-impact plastic, the SR300 Series is as durable as it is functional.

The SR300 Series is available up to 5000 gram weight capacity with 1 gram resolution (model SR320) and 500 gram weight capacity with a 0.1 gram resolution (model SR325).



2000 grams x 1 gram resolution.

**SR320** 5000 grams x 1 gram resolution.

**SR310** 200 grams x 0.1 gram resolution.

SR325 500 grams x 0.1 gram resolution.

#### **BATTERY OPERATED**

4 "AA" batteries. Auto Power Down (APD) approximately 120 seconds.

**OPTIONAL EQUIPMENT** 



Dietary portion scale



Diaper scale



SR<sup>®</sup>Instruments, Inc.

600 Young Street, Tonawanda, New York 14150 USA T: 716.693.5977 · F: 716.693.5854 · TF: 1.800.654.6360 www.srscales.com email: info-request@srinstruments.com

## Gram Scales Model SR300/SR310/SR320/SR325




### SYSTEM SPECIFICATIONS


MAX WEIGHT	<b>SR300:</b> 2000g / 70 oz <b>SR320:</b> 5000g / 175 oz <b>SR310:</b> 200g / 7 oz <b>SR325:</b> 500g / 17 oz	
DISPLAY RESOLUTION	SR300: 1g / 0.1oz SR320: 1g / 0.1oz SR310: 0.1g / 0.01oz SR325: 0.1g / 0.01oz	
DISPLAY UNITS	Grams / Ounces	
POWER	4 "AA" alkaline batteries	
AUTO POWER DOWN (APD)	Approximately 120 seconds to conserve battery life. A low battery indication will appear when it's time to replace the batteries.	
WEIGHT	2 Lb (1 Kg)	
WARRANTY	2 years	



ACCESSORIES AND OPTIONS	
DESCRIPTION	PART NUMBER
2000g SCALE	SR300
200g SCALE	SR310
5000g SCALE	SR320
500g SCALE	SR325
PLASTIC BOWL	TH200B
AC ADAPTER	FRADATH004

**DIMENSIONS** in [cm]





### PACKAGING

