

Quality. Reliability. Accuracy.

Wall Mounted Height Rod

Model SR8591

The SR8591 Wall Mounted Height Rod is the only height measurement device you will ever need. With a measurement range of 24 – 78 inches, you can easily accommodate patients whether you're in a doctor's office, clinic, or hospital. The sleek telescoping design with flush folding headpiece makes this height rod as aesthetically pleasing as it is functional.

SYSTEM	
SPECIFICATIONS	
MEASUREMENT	

RANGE	24 – 78-3/4 in / 60 – 200 cm
GRADUATIONS	1/16 in / 0.1 cm
DIMENSIONS	2-1/4 x 31 x 1-1/4 in 6 x 79 x 3 cm
WEIGHT	2 lb / 1 kg

ALUMINUM CONSTRUCTION

Durable and easy to clean.


FULL RANGE HEIGHT MEASUREMENT

Measurement range: 24 - 78 in / 60 - 200 cm.

SECURE WALL MOUNT

The height rod is secured to the wall using two low-profile rubber bumpers. All mounting hardware is included.

PACKAGING



Specifications subject to change without notice.