

The SRV415 multi-purpose mobile vet scale is a sleek, innovative scale that provides years of reliable service and accurate weight data. Designed to be low profile and durable, this scale is ideal for small to medium size animals where space is a premium, yet a weighing solution is needed. The large, comfortable platform measures 26" W x 18" L and yet can be easily picked up and moved throughout a facility. The large capacity makes it also possible for owners to hold their pets and still get accurate weight data on the animal. The scale's built-in readout makes it ideal for mobile veterinary practices.

Built-in Display

Enables scale to be moved effortlessly between locations or exam rooms.

Low Profile Design

Easy and safe access to weighing surface.

Solid Cast Aluminum Construction

Durable, stable, and easy to clean.

Contoured Base

Unique, beveled edges slope upwards to scale platform. The entire top of the scale, including the beveled edges, are part of the weighing surface.

Memory Recall

Automatically stores the last stable weight.




SYSTEM SPECIFICATIONS

Max Weight	500 Lb/221 Kg
Display Resolution	0.1 Lb / .1 Kg
Accuracy	0.1% +/-1 digit of displayed resolution
Button Functions	Zero/Weigh, Lb/Kg, Hold/Recall
Display Units	Pounds/Kilograms
Hold	Freezes the weight on the display
Recall	Recalls the last stored weight
Power	6 "AA" batteries
Portability	Built in handle
Auto Zero	One button operation
Warranty	Two (2) Years


Multi-Purpose Mobile Vet Scale

Model SRV415

DIMENSIONS in[cm]



PACKAGING



SR[®] Instruments, Inc. 600 Young Street, Tonawanda, NY 14150 USA

TF: 1.888.WEIGH.ME Ph: 716.693.5977 F: 716.693.5854

www.srscales.com email: info-request@srscales.com