

# SR Scales®

*“The Measure of Excellence”*

## Zoological Platform Scale Model **SR**®V712W

**W**hen it comes to zoo animals, SR has the leading edge technology to weighing penquins, sea otters and other small animals. Our scales provide the option for allowing animals to be weighed in their habitat or holding cage.

Constructed of heavy duty stainless steel, with sealed load cells, the SR Zoological Platform Scale provides years of dependable service.

The **SR**®V712W built in readout is equipped with bluetooth wireless capability allowing the scale to be operated remotely with the free **SR Scales**® app.

The **SR**®V712W comes with built in transport wheels and handle to allow flexibility and easy movement between exhibits.



Whatever your weighing needs may be, **SR Scales**® is “The Measure of Excellence”

**Capacity-** 1,000 pound capacity (454 Kgs)

**Heavy Duty-** Stainless Steel Platform

**Platform Size -** SRV712W 22" x 40" (56 x 102cm)  
SRV712W-2 26" x 44" (66 x 112cm)

**App-enabled -** Scale comes with Tablet preloaded with **SR Scales**® app allowing remote operation

**Battery Operated** Six "AA" Cell Batteries



Wireless App shown on Tablet



**SR**® Instruments, Inc. 600 Young Street, Tonawanda, NY 14150 USA  
TF: 1.888.WEIGH.ME Ph: 716.693.5977 F: 716.693.5854  
www.srscales.com email: info@srscales.com

MADE IN THE USA   
171012

## SYSTEM SPECIFICATIONS

Max Weight	1,000lb/ 454kg
Display Resolution	.2lb/.1kg
Display Units	Pounds or Kilograms

<u>Wireless App</u>	
Four Buttons	Zero, Weigh, Recall, and Mode

<u>Built in readout</u>	
Two button	one press weigh, hold for zero, mode

Portability	Built-in wheels and handle to provide mobility
-------------	--

Scale Weight	
SRV712W	65Lb./30kg.
SRV712W-2	80Lb./36Kg.

Power	Six "AA" Cell Batteries
-------	-------------------------


Auto Power Down	Approximately ten (10) Minutes of Non-Use
-----------------	---

Warranty	Two (2) Years
----------	---------------

# Zoological Platform Scale

Model **SR**<sup>®</sup>V712W

DIMENSIONS in (cm)



Shipping weight

