

Quality. Reliability. Accuracy.

# Table Top Scale Model SRV930

When it comes to weighing cats or smaller animals, the SRV930 is the scale for the job!

Constructed of high impact plastic, this scale is as durable as it is functional.

The state of the art load cell technology permits technicians to obtain accurate and consistent weight even if the animal will not sit perfectly still.

#### **CAPACITY**

Weighs animals up to 44 Lbs.

#### **MULTIPLE MODES**

Lbs/ lbs - Oz/ Kg/ and Oz.

#### LIGHTWEIGHT

Weighs only four (4) Lbs.

#### **PLATFORM SIZE**

24.5" x 14" (with extensions)

#### **MEMORY RECALL**

Recall the last stored weight.

#### **AUTO-POWER SHUTDOWN**

Programable.

#### **BATTERY OPERATED**

One (1) 9 volt battery.

#### **OPTIONAL AC POWER SUPPLY**

Scale will stay active when plugged in.





# Table Top Scale

## **Model SRV930**


SRSGAJQS<sup>®</sup> by SR Instruments, Inc.

Quality. Reliability. Accuracy.


# SYSTEM SPECIFICATIONS

MAX WEIGHT	44 Lbs. / 20 Kgs. / 704 ozs.
DISPLAY RESOLUTION	.02 lb / .01 kg / .5oz
BUTTON FUNCTIONS	On/Off/Tare, Unit, and Memory
DISPLAY UNITS	Pounds / ounces / Kilograms / Pounds w ounces
HOLD	Freezes the weight on the display
MEMORY	Recalls the last stored weight
TARE	Allows objects to be placed on scale and zeroed out.
POWER	One (1) 9 volt battery.
OPTIONAL POWER	AC adaptor. Scale will stay on and ready for use when plugged in.
AUTO POWER DOWN	Approximately One (1) minute of non-use.
SCALE WEIGHT	Scale weighs only 5.5 Lbs. and has integrateded handles for easy lifting.
WARRANTY	Two (2) Years

## DIMENSIONS in [cm]



### **PACKAGING**



## **ACCESSORIES AND OPTIONS**

DESCRIPTION	PART NUMBER
AC ADAPTOR	SRV930AC
SOFT SIDED CARRY CASE	SRV930CC

Specifications subject to change without notice.