

Mobile Platform Scale

Model SRV947R

The low profile SRV947R is a sleek, innovative scale that provides years of reliable service and accurate weight data.

Constructed of high strength aluminum, the SRV947R is as durable as it is functional. The scale surface has an extremely durable polyurea coating that offers grip along with ease of cleaning.

With the scales AC power adaptor the scale is always on and ready to weigh without pressing any buttons on the display. Optionally, the scale can be remotely operated using the SR Scales app on your phone or tablet.

NON-SLIP PLATFORM

Durable and easy to clean polyurea coating over stainless with aggregate mixed in to create a surface animals are comfortable walking on.

PLATFORM SIZE

24" x 48" low profile platform.

AUTO-WEIGH

The scale has an auto-weigh feature that allows the scale to wake and display a weight when an animal steps on.

OPTIONAL EQUIPMENT

Bluetooth module required for use with SR Scales App (must be factory installed).



Mobile Platform Scale

Model SRV947R


SR Scales[®]
by SR Instruments, Inc.

Quality. Reliability. Accuracy.


SYSTEM SPECIFICATIONS

MAX WEIGHT	400 lb / 182 kg
DISPLAY RESOLUTION	.2 lb / .1 kg
BUTTON FUNCTIONS	Zero/Weigh, lb/kg, Hold/Recall
DISPLAY UNITS	Pounds / Kilograms
MODEL SIZE	24" x 48"
HOLD	Freezes the weight on the display
POWER	AC power adaptor
AUTO ZERO	One button operation
SCALE WEIGHT	36 Pounds
WARRANTY	Two (2) Years

DIMENSIONS in [cm]



PACKAGING



ACCESSORIES AND OPTIONS

DESCRIPTION	PART NUMBER
REPLACEMENT AC ADAPTOR	FRGSM06U09P1J
24" X 48" WIRELESS	SRV947R-BT

Please note the wireless model still includes the wall mounted display.

Specifications subject to change without notice.